**vetproof® SL *Mycoplasma* Multiplex Detection Kit  
 (*M. gallisepticum*, *M. synoviae*)  
 -5'Nuclease- real-time PCR kit**

**Order No. V 740 01**

**Instruction For Use  
 Version 2, July 2017 - EN  
 For Research Use Only (RUO)**

The **vetproof® SL *Mycoplasma* Multiplex Detection Kit (*M. gallisepticum*, *M. synoviae*)** can be used for the qualitative detection of *Mycoplasma gallisepticum* (FAM-channel) and *Mycoplasma synoviae* (VIC/HEX-channel).

**A. Real-time PCR Time-Temperature Protocol**

The following procedure is optimized for a real-time PCR instrument with FAM (for *M. gallisepticum*), VIC/HEX (for *M. synoviae*) and Cy5 (for Internal Control) detection channel. Program the PCR instrument before preparing the reaction mixes. Use the following real-time PCR-protocol for the **vetproof® SL *Mycoplasma* Multiplex Detection Kit** (for details on how to program the experimental protocol, see the Instrument Operator's Manual of your real-time PCR-cycler):

Temperature	Time	Cycle
95 °C	10 min	1
95 °C	15 sec	40
60 °C *	60 sec	

\* Fluorescence detection in this step

**Notes:** For users of the LightCycler® 480 Systems I and II: A Color Compensation is necessary and will be supplied by BIOTECON Diagnostics for users of the LightCycler® 480 Systems I and II. Please contact BIOTECON Diagnostics for further information.

**Store at -15°C to -25°C. Expiry Date: see kit label.**

**B. Preparation of the PCR Mix**

To prevent the risk of contamination with foreign DNA, we recommend that all experiment steps be performed in a PCR clean room or separated area. Filter tips are recommended for each step. Always wear gloves when handling the PCR-vessels.

1. Thaw the kit components. Using ice or lab top cooler is recommended during experiment for maintaining the enzyme activity. Briefly centrifuge the tubes to remove air bubble and drops from the inside of the cap.
2. Total reaction volume is 20 µl, the volume of DNA sample is 6 µl. Prepare a reaction mixture according to the table below. Mix the reagents in the PCR reaction tubes by pipetting up and down.

The volumes indicated below are based on a single 20 µl standard reaction. Prepare the PCR mixture by multiplying the volume by the number of reactions required plus one or two additional reactions to cover pipetting losses.

Component	Volume
Primer / Probe Mixture	4 µl
2x Real-time PCR Master Mix	10 µl
Total	14 µl

- Pipet 14 µl PCR mix into each PCR vessel.
  - For the **samples of interest**, add 6 µl sample DNA.
  - For the **negative control**, add 6 µl H<sub>2</sub>O PCR-grade instead of sample DNA.
  - For the **positive control**, add 6 µl of Control DNA instead of sample DNA.
- Seal the PCR vessels accurately with optical caps or sealing foil.
- Briefly spin the PCR vessels in a suitable centrifuge.
- Cycle the samples as described above.

### C. Data Interpretation

The amplification of *Mycoplasma gallisepticum* DNA is analyzed in the FAM detection channel. The amplification of *Mycoplasma synoviae* DNA is analyzed in the VIC/HEX, and the Internal Control in the Cy5 detection channel.

Channel FAM	Channel VIC/HEX	Channel Cy5	Result Interpretation
Positive	Negative	Positive or Negative	Positive for <i>Mycoplasma gallisepticum</i>
Negative	Positive	Positive or Negative	Positive for <i>Mycoplasma synoviae</i>
Negative	Negative	Positive	Negative for <i>Mycoplasma gallisepticum</i> and for <i>Mycoplasma synoviae</i>
Negative	Negative	Negative	Invalid

### D. Explanation of Symbols

- REF** Product number
-  Contains sufficient for "n" tests/reactions (rxns)
-  Storage Temperature
-  Batch code
-  Expiry Date

BIOTECON Diagnostics GmbH  
 Hermannswerder 17  
 14473 Potsdam – Germany  
 Phone +49 (0) 331 2300-200  
 Fax +49 (0) 331 2300-299  
 bcd@bc-diagnostics.com  
 www.bc-diagnostics.com

vetproof® is a registered trademark of BIOTECON Diagnostics. LightCycler® is a trademark of Roche Diagnostics.

# vetproof® SL *Mycoplasma* Multiplex Detection Kit (*M. gallisepticum*, *M. synoviae*)

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